




MMP7 Antibody

Product Code	CSB-PA07455A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09237
Immunogen	Recombinant Human Matrilysin protein (95-267AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrilysin (EC 3.4.24.23) (Matrin) (Matrix metalloproteinase-7) (MMP-7) (Pump-1 protease) (Uterine metalloproteinase), MMP7, MPST1 PUMP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MMP7

Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA007455A0Rb at dilution of 1:100