





PRKCSH Antibody

Product Code	CSB-PA02925A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14314
Immunogen	Recombinant Human Glucosidase 2 subunit beta protein (15-302AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucosidase 2 subunit beta (80K-H protein) (Glucosidase II subunit beta) (Protein kinase C substrate 60.1 kDa protein heavy chain) (PKCSH), PRKCSH, G19P1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PRKCSH

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA02925A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA02925A0Rb at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.