



DRAP1 Antibody

Product Code	CSB-PA05144A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14919
Immunogen	Recombinant Human Dr1-associated corepressor protein (4-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dr1-associated corepressor (Dr1-associated protein 1) (Negative cofactor 2-alpha) (NC2-alpha), DRAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DRAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.