



CTSB Antibody, HRP conjugated

Product Code	CSB-PA06974B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07858
Immunogen	Recombinant Human Cathepsin B protein (82-333AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cathepsin B (EC 3.4.22.1) (APP secretase) (APPS) (Cathepsin B1) [Cleaved into: Cathepsin B light chain; Cathepsin B heavy chain], CTSB, CPSB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CTSB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.