



RPS25 Antibody

Product Code	CSB-PA00645A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62851
Immunogen	Recombinant Human 40S ribosomal protein S25 protein (1-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S25 (Small ribosomal subunit protein eS25), RPS25
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPS25