




RPS9 Antibody

Product Code	CSB-PA00695A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46781
Immunogen	Recombinant Human 40S ribosomal protein S9 protein (2-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S9 (Small ribosomal subunit protein uS4), RPS9
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPS9

Image



Western blot

All lanes: RPS9 antibody at 2µg/ml

Lane 1: 293T whole cell lysate

Lane 2: EC109 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 21 kDa

Observed band size: 21 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.