



SNX3 Antibody

Product Code	CSB-PA01324A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60493
Immunogen	Recombinant Human Sorting nexin-3 protein (3-162AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sorting nexin-3 (Protein SDP3), SNX3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SNX3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.