




TF Antibody

Product Code	CSB-PA00250E0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02787
Immunogen	Human Transferrin (Native Protein)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serotransferrin (Transferrin) (Beta-1 metal-binding globulin) (Siderophilin), TF
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TF

Image



Western blot

All lanes: Human Transferrin antibody at 2µg/ml
+ EC109 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/15000 dilution
Predicted band size: 78 kDa
Observed band size: 78 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.