




RPL17 Antibody

Product Code	CSB-PA02504A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18621
Immunogen	Recombinant Human 60S ribosomal protein L17 protein (2-184AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L17 (60S ribosomal protein L23) (Large ribosomal subunit protein uL22) (PD-1), RPL17
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPL17

Image



Western blot
 All lanes: RPL17 antibody at 2µg/ml + EC109 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/15000 dilution
 Predicted band size: 22, 18, 27 kDa
 Observed band size: 27, 80 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.