



ANXA5 Antibody, HRP conjugated

Product Code	CSB-PA06384B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08758
Immunogen	Recombinant Human Annexin A5 protein (2-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Annexin A5 (Anchoring CII) (Annexin V) (Annexin-5) (Calphobindin I) (CBP-I) (Endonexin II) (Lipocortin V) (Placental anticoagulant protein 4) (PP4) (Placental anticoagulant protein I) (PAP-I) (Thromboplastin inhibitor) (Vascular anticoagulant-alpha) (VAC-alpha), ANXA5, ANX5 ENX2 PP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ANXA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.