




rpsB Antibody, HRP conjugated

Product Code	CSB-PA08644B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A7V0
Immunogen	Recombinant Escherichia coli 30S ribosomal protein S2 protein (2-241AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	30S ribosomal protein S2 (Small ribosomal subunit protein uS2), rpsB
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	rpsB

Image



Western Blot

Positive WB detected in: Rosseta bacteria lysate at 40µg, 20µg, 10µg, 5µg

All lanes: rpsB antibody; HRP conjugated at 0.42µg/ml

Predicted band size: 27 kDa

Observed band size: 27 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.