



rpsB Antibody, Biotin conjugated

Product Code	CSB-PA08644D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A7V0
Immunogen	Recombinant Escherichia coli 30S ribosomal protein S2 protein (2-241AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	30S ribosomal protein S2 (Small ribosomal subunit protein uS2), rpsB
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	rpsB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.