




Hemoglobin Antibody

Product Code	CSB-PA00480E0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Guinea Pig Hemoglobin protein (Native Protein)
Raised In	Rabbit
Species Reactivity	Guinea Pig
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemoglobin
Immunogen Species	Guinea Pig
Research Area	Others

Image



Western blot
 All lanes: Guinea Pig Hemoglobin antibody at 2µg/ml
 Lane 1: Guinea Pig serum at 1: 100
 Lane 2: Guinea Pig serum at 1: 1000
 Secondary
 Goat polyclonal to rabbit IgG at 1/15000 dilution
 Predicted band size: 16 kDa
 Observed band size: 16 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.