



Csf1 Antibody

Product Code	CSB-PA07777A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07141
Immunogen	Recombinant Mouse Macrophage colony-stimulating factor 1 protein (33-262AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Macrophage colony-stimulating factor 1 (CSF-1) (MCSF) [Cleaved into: Processed macrophage colony-stimulating factor 1], Csf1, Csfm
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Csf1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.