




CYB5A Antibody

Product Code	CSB-PA03134A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00167
Immunogen	Recombinant Human Cytochrome b5 protein (2-134AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome b5 (Microsomal cytochrome b5 type A) (MCB5), CYB5A, CYB5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CYB5A

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA03134A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.