



Recombinant Human Hepatocyte nuclear factor 1-alpha (HNF1A)

Product Code	CSB-EP010592HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P20823
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVSKLSQLQT ELLAALLESGLSKEALIQALGEPGPYLLAGEPLDKGESC GGGRGELAEPLNGLGETRGS EDETDDDDGEDFTPPILKELE NLSPEEAAHQ KAVVETLLQEDPWRVAKMVK SYLQQHNIPQ REVVDTTGLN QSHLSQHLNK GTPMKTQKRA ALYTWYVRKQ REVAQQFTHA GQGGLIEEPT GDELPTKKGR RNRFKWGPAS QQILFQAYER QKNPSKEERE TLVEECNRAE CIQRGVSPSQ AQGLGSNLVTEVRVYNWFAN RRKEEAFRHK LAMDTYSGPP PGPGPGPALP AHSSPGLPPP ALSPSKVHGV RYGQPATSET AEPVSSSSGGP LVTVSTPLHQ VSPTGLEPSH SLLSTEAKLV SAAGGPLPPV STLTALHSLE QTSPGLNQQP QNLIMASLPG VMTIGPGEPA SLGPTFTNTG ASTLVIGLAS TQAQSVPVIN SMGSSLTTLQ PVQFSQPLHP SYQQPLMPPV QSHVTQSPFM ATMAQLQSPH ALYSHKPEVA QYTHTGLLPQ TMLITDTTNL SALASLTPTK QVFTSDTEAS SEGLHTPAS QATTLHVPSQ DPAGIQHLQP AHRLSASPTV SSSSLVLYQS SDSSNGQSHL LPSNHSVIET FISTQMASSS Q
Source	E.coli
Target Names	HNF1A
Protein Names	Recommended name: Hepatocyte nuclear factor 1-alpha Short name= HNF-1-alpha Short name= HNF-1A Alternative name(s): Liver-specific transcription factor LF-B1 Short name= LFB1 Transcription factor 1 Short name= TCF-1
Expression Region	1-631
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This protein is a transcription factor required for the expression of several liver-specific genes. The encoded protein functions as a homodimer and binds to the inverted palindrome 5'-GTTAATNATTAAC-3'. Defects in this gene are a cause of maturity onset diabetes of the young type 3 (MODY3) and also can result in the appearance of hepatic adenomas.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.