




IgG Antibody

Product Code	CSB-PA00770E0Gp
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Goat IgG heavy chain (Native Protein)
Raised In	Guinea pig
Species Reactivity	Goat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	IgG
Immunogen Species	Goat
Research Area	Others

Image



Western blot

All lanes: Goat IgG heavy chain antibody at 2µg/ml + Goat serum at 1: 100

Secondary

Rabbit polyclonal to Guinea pig IgG at 1/15000 dilution

Predicted band size: 55 kDa

Observed band size: 55 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.