




IgG Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA00770E0Gp |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Immunogen | Goat IgG heavy chain (Native Protein) |
| Raised In | Guinea pig |
| Species Reactivity | Goat |
| Tested Applications | ELISA, WB; Recommended dilution: WB:1:1000-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | IgG |
| Immunogen Species | Goat |
| Research Area | Others |

Image



Western blot
All lanes: Goat IgG heavy chain antibody at 2µg/ml + Goat serum at 1: 100
Secondary
Rabbit polyclonal to Guinea pig IgG at 1/15000 dilution
Predicted band size: 55 kDa
Observed band size: 55 kDa