



IgG Antibody, FITC conjugated

Product Code	CSB-PA00770G0Gp
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Goat IgG heavy chain (Native Protein)
Raised In	Guinea pig
Species Reactivity	Goat
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	IgG
Immunogen Species	Goat
Research Area	Others