




Tf Antibody

Product Code	CSB-PA00470E0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12346
Immunogen	Rat Transferrin (Native Protein)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serotransferrin (Transferrin) (Beta-1 metal-binding globulin) (Liver regeneration-related protein LRRG03) (Siderophilin), Tf
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Tf

Image



Western blot
 All lanes: Rat Transferrin antibody at 2µg/ml
 Lane 1: Rat Transferrin at 0.2 mg
 Lane 2: Rat Transferrin at 0.1 mg
 Secondary
 Goat polyclonal to rabbit IgG at 1/15000 dilution
 Predicted band size: 77, 55 kDa
 Observed band size: 77 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.