



ALB Antibody

Product Code	CSB-PA00070E0Gp
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49065
Immunogen	Serum Albumin purified from Rabbit serum (Native Protein)
Raised In	Guinea pig
Species Reactivity	Rabbit
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin, ALB
Immunogen Species	Rabbit
Research Area	Cardiovascular
Target Names	ALB