



Hemoglobin Antibody

Product Code	CSB-PA00790E0Gp
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Pigeon Hemoglobin (Native Protein)
Raised In	Guinea pig
Species Reactivity	Pigeon
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemoglobin
Immunogen Species	Pigeon
Research Area	Others