



TTR Antibody, HRP conjugated

Product Code	CSB-PA00340F0Gp
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02766
Immunogen	Human Prealbumin protein (Native Protein)
Raised In	Guinea pig
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transthyretin (ATTR) (Prealbumin) (TBPA), TTR, PALB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TTR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.