



TTR Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA00340G0Gp |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P02766 |
| Immunogen | Human Prealbumin protein (Native Protein) |
| Raised In | Guinea pig |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Transthyretin (ATTR) (Prealbumin) (TBPA), TTR, PALB |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | TTR |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |