



TADA3 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA03494A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O75528 |
| Immunogen | Recombinant Human Transcriptional adapter 3 protein (1-238AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Transcriptional adapter 3 (ADA3 homolog) (hADA3) (STAF54) (Transcriptional adapter 3-like) (ADA3-like protein), TADA3, ADA3 TADA3L |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | TADA3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |