



ENO2 Antibody, Biotin conjugated

Product Code	CSB-PA09764D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09104
Immunogen	Recombinant Human Gamma-enolase protein (2-431AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-enolase (EC 4.2.1.11) (2-phospho-D-glycerate hydro-lyase) (Enolase 2) (Neural enolase) (Neuron-specific enolase) (NSE), ENO2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ENO2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.