



# Cxcl1 Antibody

<b>Product Code</b>	CSB-PA09217A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P12850
<b>Immunogen</b>	Recombinant Mouse Growth-regulated alpha protein (29-96AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Growth-regulated alpha protein (C-X-C motif chemokine 1) (Platelet-derived growth factor-inducible protein KC) (Secretory protein N51) [Cleaved into: KC(5-72) (Hematopoietic synergistic factor) (HSF) (KC-T)], Cxcl1, Gro Gro1 Mgsa Scyb1
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Immunology
<b>Target Names</b>	Cxcl1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.