




HAGH Antibody

Product Code	CSB-PA11195A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16775
Immunogen	Recombinant Human Hydroxyacylglutathione hydrolase, mitochondrial protein (50-308AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydroxyacylglutathione hydrolase, mitochondrial (EC 3.1.2.6) (Glyoxalase II) (Glx II), HAGH, GLO2 HAGH1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HAGH

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA11195A0Rb at dilution of 1:100