



TNFSF9 Antibody, HRP conjugated

Product Code	CSB-PA08047B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41273
Immunogen	Recombinant Human Tumor necrosis factor ligand superfamily member 9 protein (52-254AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor ligand superfamily member 9 (4-1BB ligand) (4-1BBL), TNFSF9
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TNFSF9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.