



# Recombinant Rabbit Sex hormone-binding globulin (SHBG)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP021252RB  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P15196  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Oryctolagus cuniculus (Rabbit)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | TQRAQDSPA VHLINGLGQE PIQVLTFDLT RLVKASSSFE LRTWDSEGLV<br>FYGDTSPKDD WFMLGLRDGR PEIQMHNPPWA QLTVGAGPRL<br>DDGSWHQVHV KIRGDSVLE VDGKEVLRSL QVSGTLHDKP QPVMKLA VGG<br>LLFPSSLRL PLVPALDGCL RRGSWLDPQA QISASAHASR RSCDVELQPG<br>IFFPPGTHAE FSLQDIPQPQ TEPWAFSLDL ELKPSEGSGR LLALGTPEDP<br>NWLSLHLQDQ KVVLSGMEP GLDLPLAWGL PLQLKLG VST AVLSQGSKKQ<br>ALGLPPSGLG PLLNLWAQPQ GRLFLGALPG EDSSASFCLD GLWAQGQKLD<br>MDKALNRSQD IWITHSCPSSP GNGTDTSH   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | SHBG  |
| <b>Protein Names</b>     | Recommended name: Sex hormone-binding globulin Short name= SHBG<br>Alternative name(s): Sex steroid-binding protein Short name= SBP Testis-specific androgen-binding protein Short name= ABP  |
| <b>Expression Region</b> | 32-398  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Target Details</b>    | This gene encodes a steroid binding protein that was first described as a plasma protein secreted by the liver but is now thought to participate in the regulation of steroid responses. The encoded protein binds each steroid molecule as a dimer formed from identical or nearly identical monomers. The use of alternate promoters and alternatively spliced transcripts have been described. Multiple transcript variants encoding different isoforms have been found for this gene. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.