



CXCL1 Antibody, Biotin conjugated

Product Code	CSB-PA06097D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09341
Immunogen	Recombinant Human Growth-regulated alpha protein (35-107AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth-regulated alpha protein (C-X-C motif chemokine 1) (GRO-alpha(1-73)) (Melanoma growth stimulatory activity) (MGSA) (Neutrophil-activating protein 3) (NAP-3) [Cleaved into: GRO-alpha(4-73); GRO-alpha(5-73); GRO-alpha(6-73)], CXCL1, GRO GRO1 GROA MGSA SCYB1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CXCL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.