




HSPE1 Antibody

Product Code	CSB-PA10857A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61604
Immunogen	Recombinant Human 10 kDa heat shock protein, mitochondrial protein (2-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	10 kDa heat shock protein, mitochondrial (Hsp10) (10 kDa chaperonin) (Chaperonin 10) (CPN10) (Early-pregnancy factor) (EPF), HSPE1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HSPE1

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA10857A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.