



PYGB Antibody, HRP conjugated

Product Code	CSB-PA10065B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11216
Immunogen	Recombinant Human Glycogen phosphorylase, brain form protein (6-843AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycogen phosphorylase, brain form (EC 2.4.1.1), PYGB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PYGB