



TUBB3 Antibody, Biotin conjugated

Product Code	CSB-PA03874D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13509
Immunogen	Recombinant Human Tubulin beta-3 chain protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin beta-3 chain (Tubulin beta-4 chain) (Tubulin beta-III), TUBB3, TUBB4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TUBB3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.