



UCHL1 Antibody, HRP conjugated

Product Code	CSB-PA11637B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09936
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase isozyme L1 protein (1-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase isozyme L1 (UCH-L1) (EC 3.4.19.12) (EC 6.-.-.-) (Neuron cytoplasmic protein 9.5) (PGP 9.5) (PGP9.5) (Ubiquitin thioesterase L1), UCHL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	UCHL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.