



Recombinant Mouse Bisphosphoglycerate mutase (Bpgm)

Product Code	CSB-YP002781MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15327
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	SKHKLIILR HGEGQWNKEN RFCSWVDQKL NNDGLEEARN CGRQLKALNF EFDLVFTSIL NRSIHTAWLI LEELGQEWVP VESSWRLNER HYGALIGLNR EKMALNHGEE QVRLWRRSYN VTPPIEESH PYFHEIYSDR RYKVCDVPLD QLPRSESLKD VLERLLPYWK ERIAPEILKG KSILISAHGN SSRALLKHLE GISDEDIINI TLPTGVPILL ELDENLRAVG PHQFLGNQEA IQAAIKKVDD QGKVKQGKQ
Source	Yeast
Target Names	Bpgm
Protein Names	Recommended name: Bisphosphoglycerate mutase Short name= BPGM EC= 5.4.2.4 Alternative name(s): 2,3-bisphosphoglycerate mutase, erythrocyte 2,3-bisphosphoglycerate synthase EC= 3.1.3.13 EC= 5.4.2.1 BPG-dependent
Expression Region	2-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.