



Recombinant Chicken Nucleophosmin (NPM1)

Product Code	CSB-MP015996CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P16039
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MEDSAMDMS MGPLRPQTFL FGCELKAEKE YQFKVDDEEN EHQLSLRTVT LGAGAKDELH VVEAEALDYE GNPTKVVLAS LKMSVQPTVS LGGFEITPPF VLRLKCGSGP VYVSGQHLVA LEEPESEDE EEDTKIGNAS TKRPASGGGA KTPQKKPKLS EDEDDDEDE DDEDDDEDDL DDDEEEIKTP MKKPAREPAG KNMQKAKQNG KDSKPSTPAS KTKTPDSKGD KSLTPKTPKV PLSLEEIKAK MQASVDKGCS LPKLEPKFAN YVKNCFRTE QKVIQALWQW RQTL
Source	Mammalian cell
Target Names	NPM1
Protein Names	Recommended name: Nucleophosmin Short name= NPM Alternative name(s): Nucleolar phosphoprotein B23 Nucleolar protein NO38 Numatrin
Expression Region	1-294
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.