



# Recombinant *Sulfolobus islandicus* Signal recognition 54 kDa protein (srp54)

<b>Product Code</b>	CSB-MP503214FPG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	C3NHT9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Sulfolobus islandicus</i> (strain Y.N.15.51 / Yellowstone #2)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLENIRDAVR KFLTRSTPYE KAVDEFIKEL QKSLISSDVN VKLVFSLTAK IKERLNKEKP PSVLERKEWF ISIVYDELSK LFGGDKEPNV NPTKLPFIIM LVGVQGSKTTAGKLAYFY KRRGYKVGLV AADVYRPAAY DQLLQLGNQI GVPVYGEPNN QNAIEIAKKG VDTFVKNKMD IIIVDTAGRH GYGEETKLE EMKEIYEALK PDDVILVIDA SIGQKAYDLA SRFHQASPIG SIIITKMDGT AKGGGALS AV AATGATIKFI GTGEKIDELE IFNAKRYVSR ILGMDIESI LEKVKGLEEY EKIQQKMEDV MEGKGLTLR DVYAQIMALR KMGPLSKVLQ HIPGLGVMLP TPSEDQLKLG EEKIRRWLAA LNSMTYKELE NPSIIDKSRM RRIAEGSGLE VEDVRELLEW YNNMNKLLKM VKRRRGSIDK LFGGKIG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	srp54
<b>Protein Names</b>	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
<b>Expression Region</b>	1-447
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.