



# Recombinant Bovine Tyrosine 3-monooxygenase (TH)

<b>Product Code</b>	CSB-BP023470BO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P17289
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	PTPNAASPQ AKGFRRVAVSE LDAQQAEAIM SPRFVGRRQS LIQDARKERE KAEAAASSSE SAEAAAWLER DGEAVLTLF ALPPTRPPAL TRAIKVFETF EAHLHHLETR PAQPLRAGSP PLECFVRCEV PGPVVPALLS ALRRVAEDVR AAGESKVLWF PRKVSELDKC HHLVTKFDPD LDLDHPGFSD QAYRQRRKLI AEIAFQYKQG DPIPHEVYTA EETATWKEVY STLRGLYPH ACREHLEAFE LLERFCGYRE DRIPQLEDVS RFLKERTGFQ LRPAAGLLSA RDFLASLAFR VFQCTQYIRH ASSPMHSPEP ECHELLGHV PMLADRTFAQ FSQDIGLASL GVSDEEIEKL STLYWFTVEF GLCKQNGEVK AYGAGLLSSY GELLHSLSEE PEIRAFDPDA AAVQPYQDQT YQPVYFVSES FSDAKDKLRS YASRIQRPF VKFDPYTLAI DVLDSPHAIR HALDGVQDEM QALAHALNAI S
<b>Source</b>	Baculovirus
<b>Target Names</b>	TH
<b>Protein Names</b>	Recommended name: Tyrosine 3-monooxygenase EC= 1.14.16.2 Alternative name(s): Tyrosine 3-hydroxylase Short name= TH
<b>Expression Region</b>	2-491
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is involved in the conversion of tyrosine to dopamine. It is the rate-limiting enzyme in the synthesis of catecholamines, hence plays a key role in the physiology of adrenergic neurons. Mutations in this gene have been associated with autosomal recessive Segawa syndrome. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.