



# Recombinant Rabbit Lipopolysaccharide-binding protein (LBP)

<b>Product Code</b>	CSB-MP012775RB
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P17454
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TNPG LITRITDKGL EYAAREGLLA LQRKLLVTL PDSGDGFRIK HFGRAQYKFY SLKIPRFELL RGTLRPLPGQ GLSLDISDAY IHVRGSWKVR KAFLRLKNSF DLYVKGLTIS VHLVLGSESS GRPTVTTSSC SSDIQNVELD IEGDLEELN LLQSQIDARL REVLESKICR QIEEAVTAHL QPYLQTLPVT TQIDSFAGID YSLMEAPRAT AGMLDVMFKG EIFPLDHRSP VDFLAPAMNL PEAHSRMVYF SISDYVFNTA SLAYHKSGYW NFSITDAMVP ADLNIRRTTK SFRPFVPLLA NLYPNMNLLEL QGTVNSEQLV NLSTENLLEE PEMDIEALVV LPSSAREPVF RLG VATNVSA TLTLNTRKIT GFLKPGRLQV ELKESKVGGF NVELLEALLN YYILNNLYPK VNEKLAHRFP LPLLRHIQLY DLLLQTHENF LLVGANIQYR RV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	LBP
<b>Protein Names</b>	Recommended name: Lipopolysaccharide-binding protein Short name= LBP
<b>Expression Region</b>	27-482
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is involved in the acute-phase immunologic response to gram-negative bacterial infections. Gram-negative bacteria contain a glycolipid, lipopolysaccharide (LPS), on their outer cell wall. Together with bactericidal permeability-increasing protein (BPI), the encoded protein binds LPS and interacts with the CD14 receptor, probably playing a role in regulating LPS-dependent monocyte responses. Studies in mice suggest that the encoded protein is necessary for the rapid acute-phase response to LPS but not for the clearance of LPS from circulation. This protein is part of a family of structurally and functionally related proteins, including BPI, plasma cholesteryl ester transfer protein (CETP), and phospholipid transfer protein (PLTP). Finally, this gene is found on chromosome 20, immediately downstream of the BPI gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.