



Recombinant Human Insulin-like growth factor-binding protein 2 (IGFBP2)

Product Code	CSB-MP011096HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P18065
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AEVLF RCPPCTPERL AACGPPPVAP PAAVAAVAGG ARMPCAELVR EPGCGCCSVC ARLEGEACGV YTPRCGQGLR CYPHPGSELP LQALVMGEGT CEKRRDAEYG ASPEQVADNG DDHSEGGGLVE NHVDSTMNML GGGGSAGRKP LKSGMKELAV FREKVTEQHR QMGKGGKHHL GLEPKKLRP PPARTPCQQE LDQVLERIST MRLPDERGPL EHLVSLHIPN CDKHGLYNLK QCKMSLNGQR GECWCVNPNT GKLIQGAPTI RGDPECHLFY NEQQEARGVH TQRMQ
Source	Mammalian cell
Target Names	IGFBP2
Protein Names	Recommended name: Insulin-like growth factor-binding protein 2 Short name= IBP-2 Short name= IGF-binding protein 2 Short name= IGFBP-2
Expression Region	36-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.