



Recombinant Arabidopsis thaliana AT-rich interactive domain-containing protein 6 (ARID6)

Product Code	CSB-EP506609DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C0SUW7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVDTEMQEQD VTFGALVETK YLEEEPLEPE NDHNPSEIPQ PLLLGDGQAN NGHGMNGGAV GVVDSERKT RRVQMLSPKT EGENAKKRKT WLLDSEAQGT DEAGTPVEQV AFLREVEAFY KESFLEFKPP KFYGQPLNIL KLWRAVVNLG GYEVVTTNKL WRQVGESFNP PKTCTTVSYT FRNFYEKALL EYEKCLRNG ELNLPGSTLI LSSSVEKEPS SHQSGSGGRA RRDSAARAMQ GWAHQRLVGS GEVTAPAVKD KGLISTPKHK KLKSIQLQKH KQQTSMDHVV TNEADKQLAA EVVDVGPVAD WVKINVKESK DSFEIFALVP GLLRKEVRIQ SDPAGKVVIT GQPEQLDNPW GITPFKKIVD LSARIDPLHT SAVMSMHGRL FIRVPFEQ
Source	E.coli
Target Names	ARID6
Protein Names	Recommended name: AT-rich interactive domain-containing protein 6 Short name= ARID domain-containing protein 6
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.