



Recombinant Laribacter hongkongensis Orotate phosphoribosyltransferase (pyrE)

Product Code	CSB-YP506933LNO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C1D6F5
Product Type	Recombinant Protein
Immunogen Species	Laribacter hongkongensis (strain HLHK9)
Purity	>85% (SDS-PAGE)
Sequence	MSDFRQDFIR FAVEEQVLR F GEFVTKAGRP SPYFFNAGLF NHGASLLSLA RFYARSISES GIAFDMLFGP AYKGIVLAGA TAMMLAEQGR DVPFAFNRKE AKDHGEGGTL IGAPLKGRVL IIDVISAGT SVRESVEIIR ANGAEPAGVA IALDRMERGQ GELSATQEVA QKFGLPVVAI ASLDDLLGFL AGSPDLADNL TRVEAYRTQY GVR
Source	Yeast
Target Names	pyrE
Protein Names	Recommended name: Orotate phosphoribosyltransferase Short name= OPRT Short name= OPRTase EC= 2.4.2.10
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.