



Recombinant Streptococcus pneumoniae L-lactate dehydrogenase (ldh)

Product Code	CSB-YP507551FNC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C1CKW1
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pneumoniae (strain P1031)
Purity	>85% (SDS-PAGE)
Sequence	MTSTKQHKKV ILVGDGAVGS SYAFALVNQG IAQELGIIIEI PQLHEKAVGD ALDLSHALAF TSPKKIYAAQ YSDCADADLV VITAGAPQKP GETRLDLVGK NLAINKSIVT QVVESGFKGI FLVAANPVDV LTYSTWKFSG FPKERVIGSG TSLDSARFRQ ALAEKLDVDA RSVHAYIMGE HGDSEFAVWS HANIAGVNLE EFLKDTQNVQ EAELIELFEG VRDAAYTIIN KKGATYYGIA VALARITKAI LDDENAVLPL SVFQEGQYGV ENVFIGQPAV VGAHGIVRPV NIPLNDAETQ KMQASAKELQ AIIDEAWKNP EFQEASKN
Source	Yeast
Target Names	ldh
Protein Names	Recommended name: L-lactate dehydrogenase Short name= L-LDH EC= 1.1.1.27
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.