



Recombinant Escherichia coli Protein viaA (viaA)

Product Code	CSB-BP509282ENU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C4ZZ23
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12 / MC4100 / BW2952)
Purity	≥85% (SDS-PAGE)
Sequence	MLTLDTLNVM LAVSEEGEGLIE EMIIALLASP QLAVFFEKFP RLKAAITDDV PRWREALRSR LKDARVPPEL TEEVMCYQQS QLLSTPQFIV QLPQILDLLH RLNSPWAEQA RQLVDANSTI TSALHTLFLQ RWRLSLIVQA TTLNQQLLEE EREQLLSEVQ ERMTLSGQLE PILADNNTAA GRLWDMSAGQ LKRGDYQLIV KYGEFLNEQP ELKRLAEQLG RSREAKSIPR NDAQMETFRT MVREPATVPE QVDGLQQSDD ILRLLPPELA TLGITELEYE FYRRLVEKQL LTYRLHGESW REKVIERPWW HKDYDEQPRG PFIVCVDTSG SMGGFNEQCA KAFCLALMRI ALAENRRCYI MLFSTEIVRY ELSGPQGIEQ AIRFLSQQFR GGTDLASCFR AIMERLQSRE WFDADAVVIS DFIAQRLPDD VTSKVKEKLR VHQHRRFHAVA MSAHGKPGIM RIFDHIWRFD TGMRSRLLRR WRR
Source	Baculovirus
Target Names	viaA
Protein Names	Recommended name: Protein viaA Alternative name(s): VWA domain protein interacting with AAA ATPase
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.