



# Recombinant Human Neutrophil cytosol factor 2 (NCF2)

<b>Product Code</b>	CSB-MP015528HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P19878
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSLVEAISLW NEGVLAADKK DWKGALDAFS AVQDPHSRIC FNIGCMYITL KNMTEAEKAF TRSINRDKHL AVAYFQRGML YYQTEKYDLA IKDLKEALIQ LRGNQLIDYK ILGLQFKLFA CEVLYNIAFM YAKKEEWKKA EEQLALATSM KSEPRHSKID KAMECVWKQK LYEPVVIPVG KLFRPNERQV AQLAKKDYLG KATVVASVVD QDSFSGFAPL QPQAAEPPPR PKTPEIFRAL EGEAHRVLFG FVPETKEELQ VM PGNIVFVL KKGNDNWATV MFNGQKGLVP CNYLEPVELR IHPQQQPQEE SSPQSDIPAP PSSKAPGRPQ LSPGQKQKEE PKEVKLSVPM PYTLKVHYKY TVVMKTQPGL PYSQVRDMVS KKLELRLEHT KLSYRPRDSN ELVPLSEDSM KDAWGQVKNY CLTLWCENTV GDQGFDPDEPK ESEKADANNQ TTEPQLKKGS QVEALFSYEA TQPEDLEFQE GDIILVLSKV NEEWLEGECK GKVGIFPKVF VEDCATTDLE STRREV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	NCF2
<b>Protein Names</b>	Recommended name: Neutrophil cytosol factor 2 Short name= NCF-2 Alternative name(s): 67 kDa neutrophil oxidase factor NADPH oxidase activator 2 Neutrophil NADPH oxidase factor 2 p67-phox
<b>Expression Region</b>	1-526
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene encodes neutrophil cytosolic factor 2, the 67-kilodalton cytosolic subunit of the multi-protein NADPH oxidase complex found in neutrophils. This oxidase produces a burst of superoxide which is delivered to the lumen of the neutrophil phagosome. Mutations in this gene, as well as in other NADPH oxidase subunits, can result in chronic granulomatous disease, a disease that causes recurrent infections by catalase-positive organisms. Alternative splicing results in two transcript variants that encode the same protein.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.