



Recombinant *Saccharomyces cerevisiae* Long chronological lifespan protein 2 (LCL2)

Product Code	CSB-YP509785SVL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C7GQS7
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain JAY291) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	FFNFGH HQQQQQQQQQ SYEDQVLNNP CDGYLCPTDL TCVAQQKDCP CPFPKSQLKC VLPDNKFVCI SKPATHNEKF RAIYDDPVKG PKAKNKGFRD CGWVSDAYKN H
Source	Yeast
Target Names	LCL2
Protein Names	Recommended name: Long chronological lifespan protein 2
Expression Region	25-131
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.