



# Recombinant *Zygosaccharomyces rouxii* Spindle pole body component KRE28 (KRE28)

<b>Product Code</b>	CSB-MP510262ZBL
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	C5DZ48
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Zygosaccharomyces rouxii</i> (strain ATCC 2623 / CBS 732 / NBRC 1130 / NCYC 568 / NRRL Y-229) ( <i>Candida mogii</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKSDLGGEV RIKSLENQTA HYSEQALLEQ EQRVLASLRE ITQNVIAMGQ ENSLVEIKGE LESKEESELV IDPSGFQEKI DTFVELVELL KVTHLEQETL DNFLRYTISS SNLLQINSVQ DAKYVELESQ VKELEQGTLE SHKREIATK GQIKNLCQEL SMAQDSINET FLDTNALEE CDALLNELTQ LRMEKQTSEE ADTIEDDPVS QTYEDWESLQ KSKLELRLLE EETSRLQSRV ESYEDYQKRS RQLSNNDPRM LQNHKALELL VELWMTKFLP QPGISHLELF PQSRKFQFDV EPTFTVVITL ADQTTFQNVQ VYRKDAKSLV MDHGLNDEIK NSYLGTTNNIY NGLNDIIHTL QRRVQAKGSN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KRE28
<b>Protein Names</b>	Recommended name: Spindle pole body component KRE28
<b>Expression Region</b>	1-370
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.