



Recombinant *Saccharomyces cerevisiae* Spindle pole body component KRE28 (KRE28)

Product Code	CSB-MP510568SVN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	C8Z692
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MDTGSASIKD YETVLTDIED SIAVSSEEV L NNQELRLKNT LHEITSSILA INEENK FVNP LRNDESLDVE GKEVFNPKI LSAKIKEFNK LMELLKLTYL EQETLDYFFR FTLSSTKPLQ LDSEKDPQFV KLNERNVDLK EEISNVQESK IEQIKAEIQE TGHNF AEKQD LINELYLEAT GDIENCWDSL NELKNLTNKE DKNMMGEKDT ILNSSDSDDF VEETYTNWQK LLFLQKQNR LTKELKEMHE VKNQIIRKGE QSKKEDSGHL MANESELCS INLLTKFWEK HFLKGSKTT ILNFEIFTQL GKVQFEIKDM QYIIAISLSD LKRPMIKDIT ILQKAGGNIV TDIEANSKFN NKYRNNTKVQ IFEVMDDIIS ELTNE
Source	Mammalian cell
Target Names	KRE28
Protein Names	Recommended name: Spindle pole body component KRE28
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.