



# Recombinant *Saccharomyces cerevisiae* Signal transduction protein MDG1 (MDG1)

<b>Product Code</b>	CSB-EP511326SVL-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	C7GTE8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain JAY291) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<pre>MQSSLPQFTF KWPKGPEAII LTGTFDDWKG TLPMVKDPSG AFEITLPVTF DSPSSKFYFK FIVDGQWLPS KDYKVNIDEG VENNFINITEED VIKQRENGSS TLVPESAGLA VSKNAPLIEP EAEKRAKCLR KFKIKRVIKT NKQTGERSIF SQEVVELPDS EDETQQVNKT GKNADGLSGT TTIENNNGV NEEKAIKPYE ENHPKVNLVK SEG YVTDGLG KTQSSERLY ELSAEDLEKE EEEEEDEKGG GKDTSTSADA EASEDQNKEP LSKSAKFEKP EEKVPVSSIT SHAKETSVKP TGKVATETQT YETKQGAPTA AAKKIEAKKA TRPSKPKGK ETPYKGVQKN PAKNGGFFKK LAQLLK</pre>
<b>Source</b>	<i>E. coli</i>
<b>Target Names</b>	MDG1
<b>Protein Names</b>	Recommended name: Signal transduction protein MDG1 Alternative name(s): Multicopy suppressor of defective G-protein 1
<b>Expression Region</b>	1-366
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.